



# The Split and Delay Line (SDL) for MID

W. Lu<sup>1, 2</sup>, T. Noll<sup>1</sup>, T. Roth<sup>2</sup>, M. Holler<sup>3</sup>, I. Agapov<sup>2</sup>, G. Geloni<sup>2</sup>,  
J. Hallmann<sup>2</sup>, G. Ansaldi<sup>2</sup>, H. Xu<sup>2</sup>, A. Madsen<sup>2</sup> and S. Eisebitt<sup>1</sup>

<sup>1</sup>Technische Universität Berlin

<sup>2</sup>European XFEL GmbH

<sup>3</sup>Paul Scherrer Institute



Bundesministerium  
für Bildung  
und Forschung




Supported by the BMBF in the "Forschungsschwerpunkt 302: Freie Elektronen Laser".



# Outline


- Overview of the project
- Technical requirements for SDL
- Current status of SDL
  - Mechanical design
  - Crystal Cage Prototype
  - Laser Interferometer
  - In-vacuum Test-stand
  - Open Issues
- Simulation works on SDL

# Standard time structure of European XFEL




Aim of the SDL: to enable X-ray pump-probe experiments, XPCS and other fast scattering experiments with  $\Delta t < 800$  ps

# Schematic view of MID in exp. hall




# Concept of the Split and Delay Line (SDL)

The delay is achieved by detour in the upper branch



# Technical Requirements



Energy  $E$

Delay  $\Delta t$

Maximum length  $L_{\max}$

Energy bandwidth  $dE/E$

5 -10 keV

-10 to 800 ps


$\approx 2\text{m}$

$< 1\text{E-}4$




# Technical Requirements

800 ps delay trajectories at 5 keV and 10 keV for different crystals




## Geometry of SDL with Si (220) crystal





# Motion concept of upper branch



Tilt sensor




## Cable Drive Slider Advantage:

- Long distance translation
- UHV compatible
- Low heat load
- Easy to cool fixed motors

## Need to be clarified:

- Accuracy and Repeatability
- Durability of cable
- ...



Parasitic tilts are more than  $\pm 250 \mu\text{rad}$   
(measurement range)

∇  
∇


Required alignment accuracy:

0.1  $\mu\text{rad}$  in Pitch angle (vertical)


8 0.2  $\mu\text{rad}$  in Roll angle (horizontal)

Not the final decision!

# Beam Splitter




## Geometrical Splitting



Thick Si(220) crystal


## Optical Splitting



on Bragg condition


Thin Si(220) crystal  
(a few micro)

# Upper Branch crystal cage



**Fine alignment**

**Compensate parasitic motion**







# Crystal Cage Prototype

Coarse Bragg  
Alignment Cage  
 $18.8^\circ\text{-}40.2^\circ (\pm 11^\circ)$

Low cost tiny  
UHV Compatible  
Stepper Motors  
2 nm at one full step  
(with 154k:1 gear box)




Retro-reflectors for  
laser interferometer  
Work for  $\pm 11^\circ$  ?




Fine alignment stage  
Can be adjusted in  $\mu^\circ$

Laser interferometer from:  
SIOS, Renshaw, Agilent, **Zygo**, Attocube, Smaract  
The Zygo system provides linear resolution of 0.3 nm and  
8 measurement channels. The interference signal can be  
guided by fiber to the measurement electronic.






# Laser Interferometer

3-axes interferometer




Test-stand for the interferometer





## In-vacuum test for the crystal cage and laser interferometer









- 0.5 Meter Cube
- Acrylglas Glass door
- $10^{-6}$  mBar Environment
- Two linear translation stages
- One rotation stage
- Optical windows and viewpoints
- Electronic and fiber feedthroughs





- **Optical splitter, thin Si(220) crystal**  
Provider? Osaka University?
- **Channel-Cut Design**
- **X-ray Diagnostics system**  
Retractable intensity monitors  
Transparent detectors
- **Temperature stabilization system (0.1K stability)**  
Temperature of the vessel  
Evacuate heat from the motors  
Cool crystals if necessary
- **Controlling system**  
Coordinated activities ( e.g. change Bragg angle for all crystals, change delay)  
Control loop for fine adjustment of crystal cages
- ...

Input SASE pulse  
9 keV, 20 fs


Pre-mono  
Si(111) or (220)

Filtered input  
spectrum



Spectrum  
output pulsesTime structure  
output pulses



Pre-mono Si(220)  
Same photon energy  
Thin beam splitter  
Low intensity ( $\sim 1.1\%$ )  
Constant Int. ratio



Pre-mono Si(111)  
Small photon energy shift (0.6 eV)  
3x thicker beam splitter  
High intensity ( $\sim 3.5\%$ )  
Fluctuating Int. ratio


# SDL, self-seeded Simulations

Input seeded pulse  
9 keV, 20 fs




output pulses

Spectrum



Time Structure



Pre-mono Si(220)  
Same photon energy  
Thin beam splitter 1:1  
Constant Int. ratio  
High intensity ( $\sim 11.2\%$ )



*Thank you!*

---





	755-SERIES HIGH-GAIN TYPE	756-SERIES MID-RANGE TYPE
<b>TOTAL RANGE *</b>	±0.9 degree	±10 degrees
<b>RESOLUTION *</b>	0.1 microradian (0.02 arc second)	1.0 microradian (0.2 arc second)
<b>REPEATABILITY</b>	1 microradian (0.2 arc second)	2 microradians (0.4 arc second)
<b>LINEARITY</b>	1% of half span, 7% of full span (typical)	0.5% of half span, 2% of full span (typical)
<b>NATURAL FREQUENCY</b>	0.8 Hz	1.3 Hz
<b>TIME CONSTANT</b>	0.5 second	0.4 second
<b>ENVIRONMENTAL</b>	-8°C to +80°C operation and storage, -25°C version available.	-25°C to +80°C operation and storage, wide ranges available; 0 to 100% humidity
<b>TEMPERATURE COEF. (Typical)</b>	Scale Factor: +0.04%/°C, Zero: ±3 μrad/°C	Scale Factor: +0.05%/°C, Zero: ±10 μrad/°C
<b>CONNECTIONS</b>	12-inch (30 cm) wires with tinned ends or mini connectors; 1-meter Kapton wire leads on vacuum-compatible units	
<b>MATERIALS</b>	304 stainless steel enclosure. Gold anodized 6061-T6 aluminum also available on request.	
<b>SIGNAL CONDITIONING</b>	Model 781 Bench-top Unit, Model 786 Rack Mount Unit, Model 84800 Card (1-channel), Model 84828 Card (1-channel), Model 83162 Card (2-channel), or IRIS-SC Tilt Switch and Controller.	